

life-changing
technology

Oticon Government Services

Get ready for even more with Oticon More™

November 2021



NEW
STYLE



Oticon More miniRITE T

The Oticon More family is expanding with more options designed to meet your patients' needs.

- Choice of rechargeable or standard battery options available in miniRITE and CROS styles
- Choice of desktop or new portable SmartCharger
- Improved MicroShell with more retention and color options
- Increased Android™ phone compatibility
- Wireless firmware updates now available

oticon
life-changing technology

YOUR CHOICE OF CONVENIENT CHARGERS FOR OTICON MORE*



NEW

Oticon SmartCharger*

- Protective Lid
- Power Bank
- Drying Capabilities



Desktop Charger

- Easy insertion and removal**
- Contact Free**
- Quick charge**

* Oticon SmartCharger is compatible with Oticon More devices on Firmware 1.1.1
Oticon SmartCharger is not compatible with Opn S miniRITE R

**Features for both Oticon SmartCharger and Desktop charger.

OTICON MORE, OPN S, OPN, OPN PLAY & XCEED WIRELESS ACCESSORIES



ConnectClip

- Make and receive hands-free calls to both Android™ and iPhone® devices
- Wirelessly stream music, videos or other sounds from your mobile phone to your hearing aids
- Control the volume and change programs remotely
- Remote microphone



TV Adapter 3.0

- Seamlessly stream audio from a TV
- Wireless streaming connection up to 45 ft
- Digital stereo (PCM)
- Digital/RCA connection
- Can pair to multiple sets of Oticon More, Opn S™, Opn™ and Xceed devices



Remote Control 3.0

- Adjust volume
- Switch programs
- Control connectivity sources



BT D800 Bluetooth Dongle

- Small USB device to improve wireless signal on computers
- Stream video calls, movies and music from your computer
- ConnectClip required



Phone Adapter 2.0

- Turn a landline phone into a wireless Bluetooth phone
- Wireless streaming connection up to 45 ft
- ConnectClip required



EduMic

- Automatic pairing
- 3.5mm input jack
- Telecoil

RITE STYLES

Opn S miniRITE

60 85 100 105



Oticon More miniRITE R*
Opn S miniRITE R*
(lithium-ion)

60 85 100 105



Oticon More miniRITE T*
Opn S miniRITE T*
Opn Play miniRITE T*

60 85 100 105




Oticon CROS**
Oticon CROS PX**
(rechargeable)



*includes Telecoil and Toggle Switch

**Available with 60 receiver

RITE Fitting Ranges



BTE STYLES

Opn S BTE
Plus Power 13



Xceed BTE
Super Power 13



Xceed BTE
Ultra Power 675



CORDA MINIFIT THIN TUBE

*Use with Opn S BTE Plus Power 13 only




Standard 0.9 mm
Corda miniFit thin tube



Power 1.3 mm
Corda miniFit thin tube

BTE Fitting Ranges



MINIFIT RECEIVER AND EARPIECE OPTIONS



Earpieces



Now available with more options

MICROSHELL COLORS



Retention Options

Skeleton	Canal Lock
Half Skeleton	Helix Lock
Semi Skeleton	

RITE AND BTE INSTRUMENT COLORS



*Not available for Opn Play

**Available for Opn Play miniRITE T only

***Available for Oticon More miniRITE R and miniRITE T only

EARMOLD MATERIAL

- Acrylic
- Silicone
- VarioTherm® (hard material which softens from body temperature when worn)

MINIFIT EARPIECE COMPATIBILITY OVERVIEW OTICON MORE™ VS OTICON OPN S™

miniFit earpieces	Open Dome*	OpenBass Dome**	Bass Dome double vent	Bass Dome single vent*	Power Dome	Grip Tip open	Grip Tip closed
Oticon More (miniFit speaker)	-	60 85	60 85 100	-	60 85 100	60 85	60 85
Oticon Opn S, Opn and Opn Play (miniFit speaker)	60 85	-	60 85 100	60 85 100	60 85 100	60 85	60 85

*For Opn S products only **For Oticon More products only

ROES EARMOLD AND RECEIVER CODES

	ROES Receiver Code	Description	Left Part No.	Right Part No.
60 miniFit	REC2	60 miniFit Speaker, Size 0	124946	124945
	REC3	60 miniFit Speaker, Size 1	124595	124590
	REC4	60 miniFit Speaker, Size 2	124596	124591
	REC5	60 miniFit Speaker, Size 3	124597	124592
	REC6	60 miniFit Speaker, Size 4	124598	124593
	REC7	60 miniFit Speaker, Size 5	124599	124594
	85 miniFit	REC8	85 miniFit Speaker, Size 0	124968
REC9		85 miniFit Speaker, Size 1	124969	124963
REC10		85 miniFit Speaker, Size 2	124970	124964
REC11		85 miniFit Speaker, Size 3	124971	124965
REC12		85 miniFit Speaker, Size 4	124973	124966
REC13		85 miniFit Speaker, Size 5	124974	124967
100 miniFit	REC14	100 miniFit Speaker, Size 1	124585	124579
	REC15	100 miniFit Speaker, Size 2	124586	124581
	REC16	100 miniFit Speaker, Size 3	124587	124582
	REC17	100 miniFit Speaker, Size 4	124588	124583
	REC18	100 miniFit Speaker, Size 5	124589	124584


	ROES Receiver Code	Description	Left Part No.	Right Part No.
Power Mold 100	REC37	Power Mold 100 Wire, Size 1	120212	120208
	REC38	Power Mold 100 Wire, Size 2	120213	120209
	REC39	Power Mold 100 Wire, Size 3	120214	117854
	REC40	Power Mold 100 Wire, Size 4	120215	120210
	REC41	Power Mold 100 Wire, Size 5	120216	120211
	Power Mold 105	REC59	Power Mold 105 Wire, Size 1	120212
REC60		Power Mold 105 Wire, Size 2	120213	120209
REC61		Power Mold 105 Wire, Size 3	120214	117854
REC62		Power Mold 105 Wire, Size 4	120215	120210
REC63		Power Mold 105 Wire, Size 5	120216	120211

Earmold Styles & Codes	
EM9	BTE Mold
EM10	RITE Mold
EM11	Slim Tube Mold
EM12	RITE Power Mold
EM14	MicroShell



Other Features	
SnapLock 13 Tubing	Red & Blue dots
Standard BTE Tubing sizes	Soft Tip
Corda (Slim Tube) Tubing sizes	Removal String
L/R Initials	Removal Handle (post)
Red & Blue mold tips	

CUSTOM STYLES


Opn 1 Custom Styles





IIC


75 85





CIC


75 85





ITC


75 85 90 100





ITE HS

75 85 90 100



ITE FS

75 85 90 100

Custom style colors

- C005** Black*
- C004** DarkBrown
- C003** Medium Brown
- C002** Light Brown
- C001** Beige

*IIC only

*Android devices need to support ASHA to allow direct connectivity to Oticon More. Please visit oticon.com/support/compatibility for more information.



Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Questions? We're here to help!
Oticon Government Services Direct Line: **877.310.9681**