

Oticon Government Services

life-changing  
technology

OTICON | Own

# Made for your patients



## Oticon Own™ is:

**Made for a superior sound experience** featuring groundbreaking technology MoreSound Intelligence™, MoreSound Amplifier™, and MoreSound Optimizer™

**Made to be unnoticed** thanks to continuous design improvements and innovations

**Made to connect** – extensive range of connectivity options including hands-free communication from select iPhone® and iPad® devices and direct streaming from select Android™ devices\*

**Made for your patients** – five styles, five colors

\*Only available in styles with 2.4 GHz. Android devices need to support ASHA to allow direct connectivity to Own and Oticon More™. iPhone and iPad devices need to support two-way hands-free communication. Please visit [oticon.com/support/compatibility](https://www.oticon.com/support/compatibility) for more information and to see which devices are compatible.

**oticon**  
life-changing technology

## RITE STYLES

Oticon More miniRITE R\*  
(lithium-ion)

60 85 100 105



Oticon More miniRITE T (312 Zinc-air)\*  
Opn Play miniRITE T (312 Zinc-air)\*

60 85 100 105




Oticon CROS\*\*  
Oticon CROS PX\*\*  
(lithium-ion)



\*Includes Telecoil and Toggle Switch

\*\*Available with 60 or 85 receiver

### RITE Fitting Ranges



60

85

100

105

## BTE STYLES

More miniBTE R (lithium-ion)  
More miniBTE T (312 Zinc-air)

85



Opn S BTE  
Plus Power 13

105



Xceed BTE  
Super Power 13

110




Xceed BTE  
Ultra Power 675

120



### BTE Fitting Ranges



85

105

110

120

## MINIFIT RECEIVERS



60

Receiver



85

Receiver



100

Receiver



100 105

Power Receiver Mold

## CORDA MINIFIT THIN TUBE\*



Standard 0.9 mm  
Corda miniFit thin tube



Power 1.3 mm  
Corda miniFit thin tube

\*miniFit thin tubes are compatible with Oticon More miniBTE or Opn S BTE Plus Power styles only.

## MINIFIT RITE EARMOLD OPTIONS

Standard Mold

Hollow Mold

Embedded Mold



60 85

Canal



60 85

MicroMold



60 85

Hollow Canal  
LiteTip



60 85

MicroShell  
Embedded  
Receiver



100 105

Power  
Receiver Mold  
Interchangeable  
Receiver Wire

\*Retention options available with all options.

## MICROSHELL COLORS



Beige



Light Brown



Translucent  
Medium Brown



Translucent  
Dark Brown



Black



Transparent



Red



Blue


### Retention Options for RITE Earmolds

Skeleton	Canal Lock
Half Skeleton	Helix Lock
Semi Skeleton	


## RITE AND BTE INSTRUMENT COLORS



44  
Silver



92  
Steel Grey\*



91  
Silver Grey\*



90  
Chroma Beige



94  
Terracotta



93  
Chestnut Brown



63  
Black



79  
Hear Pink\*\*\*

### Opn Play only



79  
Baby Pink



57  
Power Pink



58  
Aquamarine



47  
Blue



46  
Red



48  
Emerald Green



45  
Purple

\*Not available for Opn Play

\*\*\*Available for Oticon More miniRITE R, miniRITE T, miniBTE R and miniBTE T

## MINIFIT EARPIECE COMPATIBILITY OVERVIEW OTICON MORE™ VS OTICON OPN S™


miniFit earpieces	Open Dome*	OpenBass Dome**	Bass Dome double vent	Bass Dome single vent*	Power Dome	Grip Tip open	Grip Tip closed
Oticon More (miniFit speaker)	-	60 85	60 85 100	-	60 85 100	60 85	60 85
Oticon Opn S, Opn and Opn Play (miniFit speaker)	60 85	-	60 85 100	60 85 100	60 85 100	60 85	60 85

\*For Opn S products only \*\*For Oticon More products only

**NEW**

**CUSTOM STYLES**

Oticon Own 1 Custom Styles



**Custom colors**

- C005** Black
- C004** Dark Brown
- C003** Medium Brown
- C002** Light Brown
- C001** Beige

**ROES EARMOLD AND RECEIVER CODES**

	ROES Receiver Code	Description	Left Part No.	Right Part No.
<b>60 miniFit</b>	REC2	60 miniFit Speaker, Size 0	124946	124945
	REC3	60 miniFit Speaker, Size 1	124595	124590
	REC4	60 miniFit Speaker, Size 2	124596	124591
	REC5	60 miniFit Speaker, Size 3	124597	124592
	REC6	60 miniFit Speaker, Size 4	124598	124593
	REC7	60 miniFit Speaker, Size 5	124599	124594
<b>85 miniFit</b>	REC8	85 miniFit Speaker, Size 0	124968	124961
	REC9	85 miniFit Speaker, Size 1	124969	124963
	REC10	85 miniFit Speaker, Size 2	124970	124964
	REC11	85 miniFit Speaker, Size 3	124971	124965
	REC12	85 miniFit Speaker, Size 4	124973	124966
	REC13	85 miniFit Speaker, Size 5	124974	124967
<b>100 miniFit</b>	REC14	100 miniFit Speaker, Size 1	124585	124579
	REC15	100 miniFit Speaker, Size 2	124586	124581
	REC16	100 miniFit Speaker, Size 3	124587	124582
	REC17	100 miniFit Speaker, Size 4	124588	124583
	REC18	100 miniFit Speaker, Size 5	124589	124584

	ROES Receiver Code	Description	Left Part No.	Right Part No.
<b>Power Mold 100</b>	REC37	Power Mold 100 Wire, Size 1	120212	120208
	REC38	Power Mold 100 Wire, Size 2	120213	120209
	REC39	Power Mold 100 Wire, Size 3	120214	117854
	REC40	Power Mold 100 Wire, Size 4	120215	120210
	REC41	Power Mold 100 Wire, Size 5	120216	120211
	<b>Power Mold 105</b>	REC59	Power Mold 105 Wire, Size 1	120212
REC60		Power Mold 105 Wire, Size 2	120213	120209
REC61		Power Mold 105 Wire, Size 3	120214	117854
REC62		Power Mold 105 Wire, Size 4	120215	120210
REC63		Power Mold 105 Wire, Size 5	120216	120211

Earmold Styles & Codes	
EM9	BTE Mold
EM10	RITE Mold
EM11	Slim Tube Mold
EM12	RITE Power Mold
EM14	MicroShell
EM15	Hollow/LiteTip Mold

Other Features	
SnapLock 13 Tubing	Red & Blue dots
Standard BTE Tubing sizes	Soft Tip
Corda (Thin Tube) Tubing sizes	Removal String
L/R Initials	Removal Handle (post)
Red & Blue mold tips	

**Questions? We're here to help!**  
 Oticon Government Services Direct Line: **877.310.9681**

**YOUR CHOICE OF CONVENIENT CHARGERS FOR OTICON MORE\***



Oticon miniBTE SmartCharger\*

**NEW**

**Oticon SmartChargers**

- Protective Lid
- Power Bank
- Drying Capabilities



Oticon miniRITE SmartCharger\*\*



Oticon miniBTE Charger

**Desktop Chargers**

- Easy insertion and removal
- Contact Free\*\*\*
- Quick charge\*\*\*



Oticon miniRITE Charger

\*Oticon miniBTE SmartCharger is compatible with Oticon More devices on Firmware 1.3.  
 \*\*Oticon miniRITE SmartCharger is compatible with Oticon More devices on Firmware 1.1.1 or newer.  
 Oticon miniRITE SmartCharger is not compatible with Opn S miniRITE R.  
 \*\*\*Features for both Oticon SmartCharger and Desktop charger.

**OTICON OWN\*, MORE, OPN S, OPN, OPN PLAY & XCEED WIRELESS ACCESSORIES**



**ConnectClip**

- Make and receive hands-free calls to both Android™ and iPhone® devices
- Wirelessly stream music, videos or other sounds from your mobile phone to your hearing aids
- Control the volume and change programs remotely
- Remote microphone



**TV Adapter 3.0**

- Seamlessly stream audio from a TV
- Wireless streaming connection up to 50 ft
- Digital stereo (PCM)
- Digital/RCA connection
- Can pair to multiple sets of Oticon More, Own, Opn S™, Opn™ and Xceed devices
- Can pair up to 4 TV adapters to one set of hearing aids



**Remote Control 3.0**

- Adjust volume
- Switch programs
- Control connectivity sources
- Mute



**BTD800 Bluetooth Dongle**

- Small USB device to improve wireless signal on computers
- Stream video calls, movies and music from your computer
- ConnectClip required



**Phone Adapter 2.0**

- Turn a landline phone into a wireless Bluetooth phone
- Wireless streaming connection up to 65 ft
- ConnectClip required



**EduMic**

- Automatic pairing
- 3.5mm input jack
- Telecoil
- Universal FM receiver port

\*Only available in styles with 2.4 GHz.