

Oticon Government Services

More potential than ever

life-changing
technology



NEW
STYLES



Oticon More
miniBTE R

Oticon More
miniBTE T

The Oticon More™ family is expanding with even more options designed to fit your patients' needs.

- Clinical research shows Oticon More gives the brain what it needs to make sense of sound
- NEW miniBTE styles join the existing miniRITE family to bring the superior sound quality of Oticon More to more patients
- Two-way hands free phone or video calls from select iPhone and iPad devices*

*Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad running iPadOS 15.2 or later. See which hearing aids and devices are compatible here: oticon.com/support/compatibility

oticon
life-changing technology

RITE STYLES

Opn S miniRITE Oticon More miniRITE R* (lithium-ion) Oticon More miniRITE T* Opn Play miniRITE T* Oticon CROS** Oticon CROS PX** (rechargeable)

60 85 100 105

60 85 100 105


60 85 100 105



*Includes Telecoil and Toggle Switch

**Available with 60 or 85 receiver

RITE Fitting Ranges



BTE STYLES

More miniBTE R More miniBTE T Opn S BTE Plus Power 13 Xceed BTE Super Power 13 Xceed BTE Ultra Power 675

85


105

110

120



BTE Fitting Ranges



MINIFIT RECEIVER AND EARPIECE OPTIONS



60 Receiver 85 Receiver 100 Receiver

Earpieces



60 85 MicroShell Embedded Receiver 100 105 Power Receiver Mold Interchangeable Receiver Wire MicroMold Hollow Canal LiteTip Canal Skeleton Half-Skeleton

Now available with more options

CORDA MINIFIT THIN TUBE



Standard 0.9 mm Corda miniFit thin tube

Power 1.3 mm Corda miniFit thin tube
Use with Opn S BTE Plus Power 13 only

RITE AND BTE INSTRUMENT COLORS



*Not available for Opn Play **Available for Opn Play miniRITE T only ***Available for Oticon More miniRITE R, miniRITE T, miniBTE R and miniBTE T

MICROSHELL COLORS



Beige Light Brown Translucent Medium Brown Translucent Dark Brown Black Transparent

Retention Options

| | |
|---------------|------------|
| Skeleton | Canal Lock |
| Half Skeleton | Helix Lock |
| Semi Skeleton | |

MINIFIT EARPIECE COMPATIBILITY OVERVIEW OTICON MORE™ VS OTICON OPN S™

| miniFit earpieces | Open Dome* | OpenBass Dome** | Bass Dome double vent | Bass Dome single vent* | Power Dome | Grip Tip open | Grip Tip closed |
|--|------------|-----------------|-----------------------|------------------------|-----------------|---------------|-----------------|
| Oticon More (miniFit speaker) | - | 60 85 | 60 85 100 | - | 60 85 100 | 60 85 | 60 85 |
| Oticon Opn S, Opn and Opn Play (miniFit speaker) | 60 85 | - | 60 85 100 | 60 85 100 | 60 85 100 | 60 85 | 60 85 |

*For Opn S products only **For Oticon More products only

ROES EARMOLD AND RECEIVER CODES

| | ROES Receiver Code | Description | Left Part No. | Right Part No. |
|--------------------|--------------------|-----------------------------|----------------------------|----------------|
| 60 miniFit | REC2 | 60 miniFit Speaker, Size 0 | 124946 | 124945 |
| | REC3 | 60 miniFit Speaker, Size 1 | 124595 | 124590 |
| | REC4 | 60 miniFit Speaker, Size 2 | 124596 | 124591 |
| | REC5 | 60 miniFit Speaker, Size 3 | 124597 | 124592 |
| | REC6 | 60 miniFit Speaker, Size 4 | 124598 | 124593 |
| | REC7 | 60 miniFit Speaker, Size 5 | 124599 | 124594 |
| | 85 miniFit | REC8 | 85 miniFit Speaker, Size 0 | 124968 |
| REC9 | | 85 miniFit Speaker, Size 1 | 124969 | 124963 |
| REC10 | | 85 miniFit Speaker, Size 2 | 124970 | 124964 |
| REC11 | | 85 miniFit Speaker, Size 3 | 124971 | 124965 |
| REC12 | | 85 miniFit Speaker, Size 4 | 124973 | 124966 |
| REC13 | | 85 miniFit Speaker, Size 5 | 124974 | 124967 |
| 100 miniFit | REC14 | 100 miniFit Speaker, Size 1 | 124585 | 124579 |
| | REC15 | 100 miniFit Speaker, Size 2 | 124586 | 124581 |
| | REC16 | 100 miniFit Speaker, Size 3 | 124587 | 124582 |
| | REC17 | 100 miniFit Speaker, Size 4 | 124588 | 124583 |
| | REC18 | 100 miniFit Speaker, Size 5 | 124589 | 124584 |


| | ROES Receiver Code | Description | Left Part No. | Right Part No. |
|-----------------------|--------------------|-----------------------------|---------------|----------------|
| Power Mold 100 | REC37 | Power Mold 100 Wire, Size 1 | 120212 | 120208 |
| | REC38 | Power Mold 100 Wire, Size 2 | 120213 | 120209 |
| | REC39 | Power Mold 100 Wire, Size 3 | 120214 | 117854 |
| | REC40 | Power Mold 100 Wire, Size 4 | 120215 | 120210 |
| | REC41 | Power Mold 100 Wire, Size 5 | 120216 | 120211 |
| Power Mold 105 | REC59 | Power Mold 105 Wire, Size 1 | 120212 | 120208 |
| | REC60 | Power Mold 105 Wire, Size 2 | 120213 | 120209 |
| | REC61 | Power Mold 105 Wire, Size 3 | 120214 | 117854 |
| | REC62 | Power Mold 105 Wire, Size 4 | 120215 | 120210 |
| | REC63 | Power Mold 105 Wire, Size 5 | 120216 | 120211 |

| Earmold Styles & Codes | |
|------------------------|-----------------|
| EM9 | BTE Mold |
| EM10 | RITE Mold |
| EM11 | Slim Tube Mold |
| EM12 | RITE Power Mold |
| EM14 | MicroShell |


| Other Features | |
|--------------------------------|-----------------------|
| SnapLock 13 Tubing | Red & Blue dots |
| Standard BTE Tubing sizes | Soft Tip |
| Corda (Thin Tube) Tubing sizes | Removal String |
| L/R Initials | Removal Handle (post) |
| Red & Blue mold tips | |


CUSTOM STYLES

Opn 1 Custom Styles





IIC







CIC







ITC






ITE HS





ITE FS



Custom Faceplate colors

- C005 Black*
- C004 DarkBrown
- C003 Medium Brown
- C002 Light Brown
- C001 Beige

*IIC only

*Android devices need to support ASHA to allow direct connectivity to Oticon More. Please visit oticon.com/support/compatibility for more information.



Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Questions? We're here to help!
 Oticon Government Services Direct Line: **877.310.9681**

YOUR CHOICE OF CONVENIENT CHARGERS FOR OTICON MORE*



Oticon miniBTE SmartCharger*



Oticon miniRITE SmartCharger**

NEW

Oticon SmartChargers

- Protective Lid
- Power Bank
- Drying Capabilities



Oticon miniBTE Charger



Oticon miniRITE Charger



Desktop Chargers

- Easy insertion and removal
- Contact Free***
- Quick charge

*Oticon miniBTE SmartCharger is compatible with Oticon More devices on Firmware 1.3.
 **Oticon miniRITE SmartCharger is compatible with Oticon More devices on Firmware 1.1.1 or newer.
 Oticon miniRITE SmartCharger is not compatible with Opn S miniRITE R.
 ***Features for both Oticon SmartCharger and Desktop charger.

OTICON MORE, OPN S, OPN, OPN PLAY & XCEED WIRELESS ACCESSORIES



ConnectClip

- Make and receive hands-free calls to both Android™ and iPhone® devices
- Wirelessly stream music, videos or other sounds from your mobile phone to your hearing aids
- Control the volume and change programs remotely
- Remote microphone



TV Adapter 3.0

- Seamlessly stream audio from a TV
- Wireless streaming connection up to 45 ft
- Digital stereo (PCM)
- Digital/RCA connection
- Can pair to multiple sets of Oticon More, Opn S™, Opn™ and Xceed devices



Remote Control 3.0

- Adjust volume
- Switch programs
- Control connectivity sources



BT D800 Bluetooth Dongle

- Small USB device to improve wireless signal on computers
- Stream video calls, movies and music from your computer
- ConnectClip required



Phone Adapter 2.0

- Turn a landline phone into a wireless Bluetooth phone
- Wireless streaming connection up to 45 ft
- ConnectClip required



EduMic

- Automatic pairing
- 3.5mm input jack
- Telecoil