

Your Information	Step 1: Account Information										Step 2: Patient Information									
	Ship to Account:										First Name:									
	Account Name:										Last Name:									
	Address:										Age:									
											Last 4 digits of Social Security #:									
	City:			State:			Zip:				Audiometric Information (required)									
	Phone #:										Hz:		250	500	1k	2k	3k	4k	6k	8k
	Contact Name:										Right:									
	Email (required):										Left:									
	Bill to Account:										Step 3: Fitter's Information									
Date:					Fitting Date:					Fitter's Name:										
Purchase Order #:										Fitter's Email:										

New Instrument																						
Model	Colors														Speaker							
	Diamond Black 63	Steel Grey 92	Silver Grey 91	Silver 44	Chroma Beige 90	Terracotta 94	Chestnut Brown 93	Hear Pink / Baby Pink 79	Emerald Green 48	Blue 47	Purple 45	Red 46	Aquamarine 58	Power Pink 57	Ear	Power Level		Size				
	1	2	3	4	5																	
Oticon More miniRITE T															L	60	85					
															R	100	105					
Oticon More miniRITE R (Lithium-ion)															L	60	85					
															R	100	105					
Oticon CROS PX (rechargeable)															L	60	85					
															R	60	85					
Oticon CROS															L	60	85					
															R	60	85					
Opn S1 miniRITE															L	60	85					
															R	100	105					
Opn S1 miniRITE T															L	60	85					
															R	100	105					
Opn Play 1 miniRITE T															L	60	85					
															R	100	105					

Grip Tips	L	R	Accessories	Domes										Notes		
				OpenBass Dome**		Open***		Single Vent Bass (0.8)***		Double Vent Bass (1.4)		Power				
Small, no vent	<input type="checkbox"/>	<input type="checkbox"/>	ConnectClip <input type="checkbox"/>	L	R	L	R	L	R	L	R	L	R	L	R	105 receiver is only available for use with a custom mold. Dome fitting is not an option. ** For Oticon More™ only *** For Oticon Opn S™ only
Small, 0.9mm vent	<input type="checkbox"/>	<input type="checkbox"/>	Remote Control 3.0 <input type="checkbox"/>	5 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	6 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Large, no vent	<input type="checkbox"/>	<input type="checkbox"/>	TV Adapter 3.0 <input type="checkbox"/>	8 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	10 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Large, 0.9mm vent	<input type="checkbox"/>	<input type="checkbox"/>	Phone Adapter 2.0 (Requires ConnectClip) <input type="checkbox"/>	12 mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Venting Information (required)							
Vent Size	Small	Medium	Large	Max	SAV	Pressure	No vent
L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Oticon Chargers	
Oticon More SmartCharger (not compatible with Oticon Opn S)	<input type="checkbox"/>
Oticon Desktop Charger	<input type="checkbox"/>

Special Instructions: _____

Custom RITE Earmolds	MicroShell for 60 & 85 Speakers (embedded speaker)		L	R	Hard material acrylic	Removal String	Options for MicroShell, MicroMold and LiteTip						
	Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		STD	Receiver level	<input type="checkbox"/> L	<input type="checkbox"/> 60	<input type="checkbox"/> 85			
	Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		STD	Receiver level	<input type="checkbox"/> R	<input type="checkbox"/> 60	<input type="checkbox"/> 85			
	Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			Receiver length	<input type="checkbox"/> L	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5
	Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			Receiver length	<input type="checkbox"/> R	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5
	Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			Shell and faceplate color (MicroShell only)	<input type="checkbox"/> Clear	<input type="checkbox"/> Beige	<input type="checkbox"/> Black	<input type="checkbox"/> Dark Brown	<input type="checkbox"/> Translucent	
	Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			Shell and faceplate color (MicroShell only)	<input type="checkbox"/> Med. Brown Translucent	<input type="checkbox"/> Light Brown				
	Molds for 60, 85		L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	*VarioTherm® silicone (clear only)	Removal String					
	LiteTip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		--Not Available--	<input type="checkbox"/>				STD		
	LiteTip with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		--Not Available--	<input type="checkbox"/>						
	LiteTip Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		--Not Available--	<input type="checkbox"/>						
	MicroMold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>				STD		
	Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>						
	Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>						
	Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		--Not Available--	--Not Available--						
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>							
Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		--Not Available--	--Not Available--							
Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>							
CROS Mold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>							
miniFit Power Molds for 100/105 speaker				L	R	Removal String	Options for miniFit Power Molds						
Power Receiver Mold				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Receiver level	<input type="checkbox"/> L	<input type="checkbox"/> 100	<input type="checkbox"/> 105			
Power Receiver Mold with Canal Lock				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Power Receiver wire length	<input type="checkbox"/> R	<input type="checkbox"/> 100	<input type="checkbox"/> 105			
							Power Receiver wire length	<input type="checkbox"/> L	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5
							Power Receiver wire length	<input type="checkbox"/> R	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5
							Clear shell - Faceplate color	<input type="checkbox"/> Clear	<input type="checkbox"/> Beige	<input type="checkbox"/> Light Brown	<input type="checkbox"/> Medium Brown	<input type="checkbox"/> Dark Brown	<input type="checkbox"/> Black

Please note: Scanned impressions will be kept on file for 5 years. Impression will be destroyed unless otherwise requested.

Custom Mold Styles

LiteTip	MicroMold	Power Receiver Mold	MicroShell	Skeleton	Half Skeleton	Semi Skeleton
						
(Hollow)	(Solid)		New options			(Available in hard material only)

Retention Locks

All mold styles are offered with Canal locks for better retention. Helix locks, Skeleton locks, Half Skeleton, and Semi Skeleton styles are also available upon request. (Retention Lock options vary by style and material selected.)

Canal Lock	Helix Lock
	
	(Available in hard material only)

*VarioTherm Material

VarioTherm is a body temperature-reactive plastic for soft flexible earmolds. The unique properties of this material allow the mold to remain hard at room temperature for ease of insertion and then soften (become more flexible) at body temperature for increased comfort and optimum sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer before insertion or removal of the receiver.

*miniFit 60 and 85 speakers only

Wax Protection System

Custom molds (except LiteTips) use the ProWax system. All miniFit receivers and MicroShells use the ProWax miniFit wax system.

