

Your Information	Step 1: Account Information								Step 2: Patient Information								
	Ship to Account:								First Name:								
	Account Name:								Last Name:								
	Address:								Age:								
									Last 4 digits of Social Security #:								
	City:				State:		Zip:		Audiometric Information (required)								
	Phone #:								Hz:	250	500	1k	2k	3k	4k	6k	8k
	Contact Name:								Right:								
	Email (required):								Left:								
	Bill to Account:								Step 3: Fitter's Information								
Date:				Fitting Date:				Fitter's Name:									
Purchase Order #:								Fitter's Email:									

Existing Instruments															
Please indicate existing instrument information:								Please indicate existing instrument information:							
Model: _____ Serial #: _____								Model: _____ Serial #: _____							
What type of speaker are you fitting?								What type of speaker are you fitting?							
<input type="checkbox"/> Standard <input type="checkbox"/> Medium <input type="checkbox"/> Power <input type="checkbox"/> miniFit 60 <input type="checkbox"/> miniFit 85 <input type="checkbox"/> miniFit 100 <input type="checkbox"/> miniFit 105 <input type="checkbox"/> DesignRITE 80								<input type="checkbox"/> Standard <input type="checkbox"/> Medium <input type="checkbox"/> Power <input type="checkbox"/> miniFit 60 <input type="checkbox"/> miniFit 85 <input type="checkbox"/> miniFit 100 <input type="checkbox"/> miniFit 105 <input type="checkbox"/> DesignRITE 80							

BTE Mold Order											
Type	L	R	Style	L	R	Options	L	R	Tubing	L	R
BTE Mold	<input type="checkbox"/>	<input type="checkbox"/>	Canal	<input type="checkbox"/>	<input type="checkbox"/>	Name/initials on mold:	<input type="checkbox"/>	<input type="checkbox"/>	Snap Lock tube (13 Medium)	<input type="checkbox"/>	<input type="checkbox"/>
Slim Tube Mold (Corda)	<input type="checkbox"/>	<input type="checkbox"/>	Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	_____	<input type="checkbox"/>	<input type="checkbox"/>	13 Medium tube	<input type="checkbox"/>	<input type="checkbox"/>
Material	L	R	Half Shell	<input type="checkbox"/>	<input type="checkbox"/>	Print L&R on mold	<input type="checkbox"/>	<input type="checkbox"/>	13 Heavy tube	<input type="checkbox"/>	<input type="checkbox"/>
Hard	<input type="checkbox"/>	<input type="checkbox"/>	Full Shell	<input type="checkbox"/>	<input type="checkbox"/>	Make Canal Tip red & blue	<input type="checkbox"/>	<input type="checkbox"/>	13 Double-wall tube	<input type="checkbox"/>	<input type="checkbox"/>
VarioTherm® (clear only)	<input type="checkbox"/>	<input type="checkbox"/>	Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	Blue/red/white dot on mold	<input type="checkbox"/>	<input type="checkbox"/>	Dry Tube (13 Medium)	<input type="checkbox"/>	<input type="checkbox"/>
Soft silicone (clear only)	<input type="checkbox"/>	<input type="checkbox"/>	Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	Removal String	<input type="checkbox"/>	<input type="checkbox"/>	Slim Tube Only available on Opn S 1 Plus Power	L	R
Color	L	R	Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	Removal String heavy	<input type="checkbox"/>	<input type="checkbox"/>	Corda2 miniFit tube 1.3, size -1	<input type="checkbox"/>	<input type="checkbox"/>
Clear (STD)	<input type="checkbox"/>	<input type="checkbox"/>	Length	L	R	Removal handle	<input type="checkbox"/>	<input type="checkbox"/>	Corda2 miniFit tube 1.3, size 0	<input type="checkbox"/>	<input type="checkbox"/>
Beige	<input type="checkbox"/>	<input type="checkbox"/>	Short Canal	<input type="checkbox"/>	<input type="checkbox"/>	Retention Ring	<input type="checkbox"/>	<input type="checkbox"/>	Corda2 miniFit tube 1.3, size 1	<input type="checkbox"/>	<input type="checkbox"/>
Light brown	<input type="checkbox"/>	<input type="checkbox"/>	Medium Canal	<input type="checkbox"/>	<input type="checkbox"/>				Corda2 miniFit tube 1.3, size 2	<input type="checkbox"/>	<input type="checkbox"/>
Medium brown	<input type="checkbox"/>	<input type="checkbox"/>	Long Canal	<input type="checkbox"/>	<input type="checkbox"/>				Corda2 miniFit tube 1.3, size 3	<input type="checkbox"/>	<input type="checkbox"/>
Dark brown	<input type="checkbox"/>	<input type="checkbox"/>	As Marked	<input type="checkbox"/>	<input type="checkbox"/>				Corda2 miniFit tube 1.3, size 4	<input type="checkbox"/>	<input type="checkbox"/>

Venting Information (required)								Notes	
Vent Size	Small	Medium	Large	Max	SAV	Pressure	No vent	Power Molds default to 1.4mm vents unless otherwise requested.	
L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	LiteTip and MicroMolds come with largest vent possible unless otherwise requested.	
R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Special Instructions:

Custom RITE Earmolds	MicroShell for 60 & 85 Speakers (embedded speaker)		L	R	Hard material acrylic	Removal String	Options					
	Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	STD	MicroShell (acrylic only)	<input type="checkbox"/> L	<input type="checkbox"/> 60	<input type="checkbox"/> 85			
	Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	STD		<input type="checkbox"/> R	<input type="checkbox"/> 60	<input type="checkbox"/> 85			
	Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Receiver Length	<input type="checkbox"/> L	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5
	Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/> R	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5
	Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Shell and Faceplate color	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Black <input type="checkbox"/> Dark Brown Translucent <input type="checkbox"/> Med. Brown Translucent <input type="checkbox"/> Light Brown					
	Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>								
	Molds for 60, 85	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	*VarioTherm® silicone (clear only)	Removal String					
	LiteTip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	--Not Available--	<input type="checkbox"/>	STD					
	LiteTip with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	--Not Available--	<input type="checkbox"/>						
LiteTip Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	--Not Available--	<input type="checkbox"/>							
MicroMold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	STD						
Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	--Not Available--	--Not Available--							
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	--Not Available--	--Not Available--							
Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
CROS Mold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
miniFit Power Molds for 100/105 speaker	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	*VarioTherm® silicone (clear only)	Removal String						
Power Receiver Mold	<input type="checkbox"/>	<input type="checkbox"/>	STD	--Not Available--	--Not Available--							
Power Receiver Mold with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	--Not Available--	--Not Available--							
Clear shell - Face plate color:	<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Light Brown <input type="checkbox"/> Medium Brown <input type="checkbox"/> Dark Brown <input type="checkbox"/> Black											
Power Receiver Wire Length	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5											

Custom Mold Styles

LiteTip	MicroMold	Power Receiver Mold	MicroShell	Skeleton	Half Skeleton	Semi Skeleton
						
(Hollow)	(Solid)					(Available in hard material only)

Retention Locks

All mold styles are offered with Canal locks for better retention. Helix locks, Skeleton locks, Half Skeleton, and Semi Skeleton styles are also available upon request. (Retention Lock options vary by style and material selected.)

Canal Lock	Helix Lock
	 (Available in hard material only)

*VarioTherm Material

VarioTherm is a body temperature-reactive plastic for soft flexible earmolds. The unique properties of this material allow the mold to remain hard at room temperature for ease of insertion and then soften (become more flexible) at body temperature for increased comfort and optimum sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer before insertion or removal of the receiver.

*miniFit 60 and 85 speakers only

Wax Protection System

Custom molds (except LiteTips) use the ProWax system. All miniFit receivers and MicroShells use the ProWax miniFit wax system.



21-137128 15500-2218/09.21