

# miniRITE Earmold options

## RITE MOLDS

Standard Mold	Hollow Mold	Embedded Mold	
---------------	-------------	---------------	--



**Canal**



**LiteTip**



**MicroShell**



**Power Receiver Mold**

**Power Options**

60 dB	•	•	•
85 dB	•	•	•
100 dB			•
105 dB			•

**Material**

Acrylic	•	•	•
Silicone	•		

**Style/Retention**

MicroMold (CIC)	•	•	•
Canal	•	•	•
Canal Lock	•	•	•
Half-Skeleton	•	•	•
Semi-Skeleton	•	•	•
3/4 Skeleton	•	•	•
Skeleton	•	•	•
Helix Lock	•	•	•
Half Shell	•	•	•
3/4 Shell	•	•	•
Full Shell	•	•	•

## Now more color options for MicroShell



Beige



Light Brown



Translucent  
Medium Brown



Translucent  
Dark Brown



Black



Transparent










Red



Blue

# miniRITE Instant Fit Options

When fittings do not require a custom earmold, Oticon offers a selection of instant fit options.

miniFit Earpiece Compatibility Overview Oticon More™ vs Oticon Opn S™							
miniFit earpieces							
	Open Dome	OpenBass Dome	Bass Dome double vent	Bass Dome single vent	Power Dome	Grip Tip open	Grip Tip closed
Oticon More (miniFit speaker)	-	60 85	60 85 100	-	60 85 100	60 85	60 85
Oticon Opn S, Opn and Opn Play (miniFit speaker)	60 85	-	60 85 100	60 85 100	60 85 100	60 85	60 85

Instant Fit Options							
miniFit earpieces	Open Dome*	OpenBass Dome**	Bass Dome double vent	Bass Dome single vent*	Power Dome	Grip Tip open	Grip Tip closed
Domes & Grip Tips	5mm† 6mm 8mm 10mm -	5mm† 6mm 8mm 10mm 12mm	- 6mm 8mm 10mm 12mm	- 6mm 8mm 10mm 12mm	- 6mm 8mm 10mm 12mm	Large/Small Vent	Large/Small No Vent
Grip Tips	<b>Grip Tip</b> provides an instant fit option that is more durable than domes. Soft flexible tips shape comfortably to most ears; the slightly tacky texture helps prevent slippage.						

\*For Opn S products only    \*\*For Oticon More products only    †For size 60 speaker only. Compatible with both Opn S and Oticon More devices.

20-42650 15555-10220/08.22