

Phone compatibility

The information below includes a list of ASHA compatible Android™ phones, allowing high-quality, one-way direct streaming for phones with Android 10 or later and Bluetooth® 5.0.

Please note that the following list of devices have been tested by Oticon and are compatible with Polaris™ platform instruments and later. This list is subject to change as Android operating systems are updated and new phones are released to the market* (this list is current as of October 2023).

Google	Google Pixel 3	Google Pixel 3 XL	Google Pixel 3a	Google Pixel 3a XL	Google Pixel 4
	Google Pixel 4a	Google Pixel 4a 5G	Google Pixel 5	Google Pixel 6	Google Pixel 6a
	Google Pixel 6 Pro	Google Pixel 7 Pro	Google Pixel 7		
OnePlus	OnePlus 10 Pro 5G	OnePlus 10R	OnePlus 7T	OnePlus 7T Pro	OnePlus 8T
	OnePlus 9	OnePlus Nord N20 5G	OnePlus 10T	OnePlus Nord 2T	OnePlus 9 Pro
Samsung	Samsung Galaxy A21s	Samsung Galaxy A50	Samsung Galaxy A51	Samsung Galaxy A52	Samsung Galaxy A52 5G
	Samsung Galaxy A52s 5G	Samsung Galaxy A70	Samsung Galaxy A71	Samsung Galaxy A71 5G	Samsung Galaxy F23 5G
	Samsung Galaxy M31	Samsung Galaxy Note10	Samsung Galaxy Note10+	Samsung Galaxy Note20 5G	Samsung Galaxy Note9
	Samsung Galaxy Note20 Ultra 5G	Samsung Galaxy S10	Samsung Galaxy S10+	Samsung Galaxy S10e	Samsung Galaxy S20
	Samsung Galaxy S20 5G	Samsung Galaxy S20 FE 5G	Samsung Galaxy S20 Ultra 5G	Samsung Galaxy S20+ 5G	Samsung Galaxy S21
	Samsung Galaxy S21 5G	Samsung Galaxy S21 FE	Samsung Galaxy S21 Ultra 5G	Samsung Galaxy S21+ 5G	Samsung Galaxy S22
	Samsung Galaxy S22 5G	Samsung Galaxy S22 Ultra	Samsung Galaxy S22+	Samsung Galaxy S9	Samsung Galaxy Tab S6
	Samsung Galaxy Tab S8 Ultra	Samsung Galaxy Z Flip3	Samsung Galaxy Z Fold3	Samsung Galaxy A53 5G	Samsung Galaxy Z Fold5
	Samsung Galaxy A73 5G	Samsung Galaxy A72	Samsung Galaxy A33 5G	Samsung Galaxy M52 5G	Samsung Galaxy A12 Nacho
	Samsung Galaxy S9+	Samsung Galaxy S23	Samsung Galaxy S23+	Samsung Galaxy S23 Ultra	Samsung Galaxy A23
	Samsung Galaxy A23 5G				
Other	LG V60 ThinQ	Motorola Edge 30 Pro	Motorola Edge 30 Ultra	Oppo A74 5G	Oppo A77s
	Vivo IQOO Z6 Pro	Vivo X70 Pro+	Vivo X80	Vivo X80 Pro	Vivo V25

*Please note where variants exist US market phones were tested.

ConnectClip

Phones that do not appear on this list can still stream with Oticon hearing aids by using the Oticon ConnectClip for hands-free phone calls and audio.



Apple® Device Phone Compatibility

Apple devices running iOS® 11 or later can stream phone calls and audio via one-way direct streaming with all Oticon instruments that have 2.4GHz Bluetooth®.

Hands-free communication is available for Polaris platform instruments** and later with iPhone® 11 or later, iPad Pro® 12.9-inch (5th generation and 4th generation), iPad Pro 11-inch (3rd generation and 2nd generation), iPad Air® (4th generation), and iPad mini® (6th generation), with software updates iOS 15.2 and iPadOS® 15.2 or later.

**Oticon More™ hearing aids require firmware 1.3.0 or later to enable hands-free communication functionality.

23-529923 15500-0218/1023