

Oticon Opn Play™ classroom connectivity

A full range of easy-to-use solutions for the classroom

ConnectClip remote microphone

ConnectClip connects to integrated 2.4 GHz receiver



Universal ear-level options

- Connect to Oticon Amigo FM
- Connect to systems from other manufacturers



FM10 adapter



Amigo R2 receiver



Universal neckloop receiver

Amigo Arc



Dedicated FM receiver

Amigo R12G2

Oticon Amigo FM

Manufacturer:
Oticon A/S
Kongebakken 9
DK-2765 Smorum
Denmark
www.oticon.global

Local Manufacturer & Distributor:
Audmet Canada Ltd
1600-4950 Yonge St
Toronto, ON M2N 6K1
www.oticon.ca

Oticon Opn Play™ classroom connectivity



BTE PP

miniRITE R

miniRITE T

miniRITE

Oticon Amigo FM - ear-level

Setting the system up requires FM transmitter and an FM receiver. Connect FM receiver to the hearing aid either directly or using an adapter and switch on the FM transmitter to activate the FM system.

- To attach either Oticon Amigo R12G2 or the FM 10 adapter to Oticon Opn Play BTE PP, the standard battery drawer needs to be replaced with the FM adapter battery drawer.
- The FM adapter battery drawer and Oticon Amigo R12G2 are available to match instrument colours.
- The FM adapter battery drawer can be made tamper resistant by inserting the special TAR left/right colour markers.



FM adapter battery drawer



Amigo T30/T31 FM transmitter



Amigo R12G2 Dedicated FM receiver



Amigo R2 universal FM receiver



FM 10 Universal FM adapter

Let children connect easily to the world around them

Connectivity options

| | | | | |
|----------------------------|---|---|---|---|
| Stereo streaming (2.4 GHz) | • | • | • | • |
| Oticon ON App | • | • | • | • |
| ConnectClip | • | • | • | • |
| Remote Control 3.0 | • | • | • | • |
| TV Adapter 3.0 | • | • | • | • |
| Phone Adapter 2.0 | • | • | • | • |

A full range of easy-to-use solutions for the classroom

Classroom receiver options

| | | | | |
|---------------------------------|-------------|-----------|-----------|---|
| Integrated 2.4 GHz receiver | • | • | • | • |
| Dedicated FM receiver | Amigo R12G2 | - | - | - |
| Universal ear-level FM receiver | Amigo R2 | - | - | - |
| Universal neckloop receiver | Amigo Arc | Amigo Arc | Amigo Arc | - |

Integrated 2.4 GHz solution

Completely integrated solution with the built-in 2.4 GHz Bluetooth™ low energy receiver in every Oticon Opn Play hearing aid, streams speech seamlessly from the remote microphone in **ConnectClip**.

Dedicated FM solution

The dedicated **Oticon Amigo R12G2** FM receiver integrates smoothly with the design of the Oticon Opn Play BTE PP. This design integrated solution allows connection to Oticon Amigo FM transmitters.

Universal ear-level solution

Oticon Opn Play BTE PP ensures full compatibility with existing classroom solutions through the **FM 10 adapter**. The student can use a universal ear-level solution by attaching any universal receiver to the FM 10 adapter.

Universal neckloop solution

Oticon Opn Play models with telecoil ensure compatibility with Oticon Amigo FM transmitters through the **Oticon Amigo Arc** FM neckloop receiver. This universal neckloop solution offers a common FM solution, and only requires one FM receiver per student.