

Oticon ♦ Tego

Easier hearing for children



child  
friendly  
hearing  
care

oticon  
PEOPLE FIRST

## Tego – the popular hearing solution for children

Like any parent, you want your child to grow up being happy, capable and independent. Since this journey demands the ability to hear and understand the spoken word, your child needs an advanced hearing solution.

Tego uses sophisticated sound processing systems to provide maximum access to speech. It prioritizes voices, controls unwanted noise, and keeps the sound quality clear and comfortable.

Thousands of children have already discovered the benefits of Tego – and your child can benefit too. With Tego they will have every chance to grow, develop and achieve their full potential.



Tego is designed to enhance all the nuances of sound that are important for your child, and it uses advanced digital sound processing to do this.

Tego offers:

### Speech priority

- > Tego uses the built-in DecisionMaker to automatically process sounds correctly in every situation.
- > Tego automatically finds the best way to make voices audible to your child.
- > Tego reduces noise without removing important background sounds.

### Child-friendly design

- > Tego comes in all sorts of fun colours that children will find attractive.
- > The case is lightweight, practical and comfortable to wear.
- > The battery door can be made tamper-resistant to prevent your child getting access to the battery.

### A reliable companion

- > Tego is robust and reliable, making it a perfect companion for active children.
- > High quality components – including a moisture-resistant shell - make Tego extremely resilient to wear and tear.

### Less feedback

- > Tego has an automatic feedback cancellation system, so your child can chew, change clothes and give people hugs, without any irritating whistling.
- > With less whistling your child can relax – and so can people in the immediate vicinity.

### Easy connection to FM

- > With optional FM equipment, Tego can help your child to hear even more clearly in class.
- > At home, personal FM solution can make it easier to communicate and interact.

Open here



# A truly child-friendly instrument



## Speech priority

Children's needs change during the course of the day, and Tego is built to help them cope. The automatic DecisionMaker in Tego always finds the best settings for each particular listening situation. In a restaurant, for instance, DecisionMaker will activate directional microphones to make the parent's voice stand out more clearly.

Tego also uses noise management to reduce noise and enhance speech. But it never removes background noise completely, because children need to be able to hear environmental sounds in order to remain safe and secure.

## Child-friendly design

Tego's slender, lightweight shell combines durability with functionality. The design is simple and practical, which is just what children need. Tego is comfortable to wear, easy to operate, and easy to keep clean.

## A reliable companion

Tego is robust and reliable, making it perfect for active little individuals. High quality exterior and interior components coupled with fully automatic operation ensure that Tego continues to support your child no matter how rough-and-tumble things get.

Tego's shell is made of robust, high quality materials. Inside the instrument the components are well protected to avoid moisture and temperature-related breakdowns. So your child can count on stable performance, every single day.

## Less feedback

Children with hearing loss can sometimes experience sudden whistling sounds known as feedback which can be irritating. Tego uses Dynamic Feedback Cancellation to prevent this. Whenever it detects feedback, a counter signal eliminates it before it can be heard.

With this kind of protection, activities like changing clothes or hugging people become far more enjoyable. With no sudden surprises your child can relax – and family and classmates can too.

## Easy connection to FM

Wireless FM systems offer valuable support in noisy classrooms. The teacher wears a microphone, and their voice is transmitted directly to the small FM receiver attached to the Tego hearing aid. This makes it seem as though they are standing just inches apart.

Tego works with all standard FM systems including the Oticon Amigo FM family. As a parent, you can wear the Amigo transmitter around your neck, or put it on your kitchen table to shorten the listening distance between you and your child. You can also plug Amigo into the TV. Your Hearing Care Professional can provide more information about Amigo FM products in general and about Direct Audio Input (DAI) options.

## Tamper resistant features

Tego's child-friendly features also include a tamper-resistant battery door and a volume control that can be disabled for younger children. Older children will find it easy to adjust the digital volume control and switch between listening programs.



## Cool, fun colours for children

Tego comes in all sorts of fun colours designed to appeal to children. These include Oticon's classic colours (including baby blue and baby pink) and cool colours. With such a wide choice it won't be hard to find one your child will be happy to wear.



- ◆ **Speech priority**
- ◆ **Child-friendly design**
- ◆ **A reliable companion**
- ◆ **Less feedback**
- ◆ **Easy connection to FM**

*People first*



We believe that it takes more than technology and audiology to create the best hearing instruments. That's why we put the individual needs and wishes of people with hearing loss first in our development of new hearing care solutions.