

IP30 clip

#8537295

The Gold
Standard

of Audiometric Transducers

Features

- Replacement velcro clip for insert earphones
- Robust, durable design
- Compatible with a variety of insert earphones, including IP30, E-A-RTone 3A, 5A, Etymotic ER-3A and ER-3C.



Technical Specifications

Material:

Metal clip with velcro disc

Weight complete:

approximately 5 grams

Dimensions:

W: 2 cm, H: 3 cm

Related accessories

Part number	Description
8507410	IP30 Fastener kit, Incl. velcro neckstrap and clip

D-0137277