


SONIC | radiant.



miniRITE T R
Rechargeable

miniRITE T
Non-rechargeable



Bluetooth is a registered trademark owned by Bluetooth SIG, Inc., USA.
Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.
App Store is a service mark of Apple Inc.
Android and Google Play are trademarks of Google LLC.

www.sonici.com
© 2021 Sonic Innovations, Inc. All rights reserved.
2021-07-20 | 236622 | UK | subject to change



radiant™
Clearly Brilliant

www.sonici.com



Powering the Evolution of the Personalized Hearing Experience.

Bright, crisp sound. Impeccable speech clarity. And a truly personalized sound experience. That's what you get with Radiant on the Extend Technology Platform. With more customization options than ever before, lithium-ion rechargeability and the ability to hear sound directly from your iPhone® and compatible Android™ devices, Radiant helps you make the most of every moment.



EXTEND Platform

For the Moments that Matter

Radiant*

The first hearing aid built on the **Extend Technology Platform** now includes more options to make it your own! Radiant provides a groundbreaking sound experience. You'll appreciate its high-performance systems that keep speech clear and reduce background noise, while virtually eliminating feedback. Enjoy streaming phone calls, music and other audio directly to Radiant from a variety of smartphones and tablets, or use its built-in telecoil in theaters and public places.

SmartMusic Pro

Listening to music is more than just a fun and relaxing hobby — it supports your overall health and well-being too. Tune in to your favorite songs with **SmartMusic Pro**. This unique program improves the sound quality of recorded, streamed and live music, so you can enjoy a better listening experience with your hearing aids.

Personalization

Get a better fit, right from the start. Your hearing care professional can create a personalized hearing experience that works for you.

Rechargeability

Radiant miniRITE T R hearing aids are equipped with built-in rechargeable **lithium-ion batteries**, which provide a full day's charge, even while streaming. Never fumble with another tiny battery again! Just set your hearing aids in the charger at the end of your day for a full recharge in three hours. For your on-the-go lifestyle, the **Charger Plus** provides convenient charging when away from home, with its built-in power bank.



miniRITE T R in Charger Plus

Radiant Hearing Aids. Powered by Innovation.

Noises can come from any direction, at any time, creating a distracting listening experience for hearing aid users. That's why the Radiant family of hearing aids has incorporated Sonic's advanced technologies that respond quickly to find, filter and fine-tune sounds arising from anywhere in your surroundings. Intelligent, integrated systems work together so you can hear more easily.



Find

Adaptive directional microphones pinpoint speech to keep conversations in focus.



Filter

Automatic noise reduction technologies filter out background sounds so you can enjoy a more comfortable listening experience.



Fine-Tune

A high-resolution amplification system enhances speech and minimizes ongoing distractions, making it easier for you to hear.

*Not all features and options are available on all technology levels or models. Ask your Hearing Care Provider for details.

Radiant Style Options

Radiant is available in both rechargeable (miniRITE T R) and non-rechargeable (miniRITE T) styles. Ask your hearing care professional which is best for you.



Get more control with the SoundLink 2 app



- Stream audio directly from all your compatible devices, including iPhones, iPads and Android phones.**
- Monitor battery charge, change your listening program, adjust volume, and more!



**Radiant is a Made for iPhone®, iPad®, iPod® hearing aid, compatible with devices running iOS 13 or later. Direct audio streaming for Android devices requires Android 10 or later, Bluetooth® 5.0 and an implementation of Audio Streaming for Hearing Aids (ASHA) on the Android device. For information on compatibility, please visit www.sonici.com/compatibility.

Convenient Wireless Accessories



SoundClip-A

Listen in stereo from Bluetooth-compatible devices or use as a remote microphone, remote control or wireless headset for mobile phones.



RC-A Remote Control

Adjust volume or programs with a single click.



TV-A Adapter

Transmit TV audio directly to your Radiant hearing aids.



Transmitter Receiver

CROS/BiCROS

For Single Sided Deafness, the transmitter picks up sound from the ear with poorer hearing and wirelessly transmits it to a Radiant receiver worn on the ear with better hearing.

For more information about Radiant, ask your Hearing Care Provider or visit www.sonici.com.